

## Japanese cedar pollen

### 1 . Product

|                       |  |
|-----------------------|--|
| Appearance :          | yellow powder  |
| Origin :              | <i>Cryptomeria japonica</i>  |
| Collecting method :   | suction  |
| Handling method :     | Foreign substances were removed by a sieve with a mesh size of 53 $\mu\text{m}$ .<br>Store at -30 °C |
| Place of collecting : |  |
| Production date :     |  |
| Lot :                 |  |
| Volume :              | 5 g  |

### 2 . Specification

The following requirement shall apply.

| Item              | Standard                            |
|-------------------|-------------------------------------|
| visual inspection | lemon yellow ~ ocher colored powder |

### 3 . Storage

$\leq$  - 20 °C

### 4 . Manufacture & Retail

Biostir Inc.  ピオスタ

Annex Shin-Osaka 4A, 2-9-1 Miyahara, Yodogawa-ku, Osaka  
532-0003, Japan  
TEL. 81-6-6151-3950 / FAX.81-6-6151-3940  
E-mail: info@biostir.com  
URL: <http://www.biostir.com>