

Atopic dermatitis challenging ointment

Biostir AD

Derive : Ointment including components of House Dust Mite, *Dermatophagoides farinae*

Lot No :

Character : Brown ointment

Quantity : 15 g

Storage&Handling : The ointment should be stored below -20°C . Repeated freezing and thawing should be avoided. The ointment should be dispensed and thawed at room temperature before use.

Use : Biostir AD is able to induce atopic dermatitis-like skin lesion in NC/Nga mice.
Please contact Biostir Inc. if you need “Protocols” and ask questions.

Note : Research use only.

Reference : Mina Yamamoto, Takayo Haruna, Kiyoshi Yasui, Hisashi Takahashi, Miho Iduhara, Shigeki Takaki, Masashi Deguchi and Akinori Arimura ; A Novel Atopic Dermatitis Model Induced by Topical Application with *Dermatophagoides Farinae* Extract in NC/Nga Mice. *Allergology International* 2007; 56: 139-148



Biostir Biostir, Inc.

Annex Shinosaka 4A, 2-9-1 Miyahara, Yodogawa-ku, Osaka, Osaka, Japan
tel +81-6-6151-3950 fax +81-6-6151-3940

e-mail: info@biostir.com

URL: <https://www.biostir.com/>