

Japanese cedar pollen extract (cpe002)

1. Product

Appearance :	yellow powder
Physical state :	powder, freeze-dried
Preparation :	100 g cedar pollen was extracted with 2 L of 0.125 M ammonium bicarbonate
Lot :	
Volume :	50 mg

2. Storage

store at -20°C

3. Notes

For research use only. Not for any clinical, therapeutic, or diagnostic use.

4. Manufacture & Retail

Biostir Inc.



Annex Shinosaka 4A, 2-9-1 Miyahara, Yodogawa-ku, Osaka, Osaka, Japan

tel +81-6-6151-3950 fax +81-6-6151-3940

e-mail: info@biostir.com

URL: <https://www.biostir.com/>